

## TARGEFFECT-HUVEC Transfection Protocols

**Transfections per kit: 600 transfections (24 well), 300 transfections (12 well), 150 transfections (6-well)**

Efficient gene transfer into primary endothelial cells is an important prerequisite for studying regulation of gene expression in vascular tissue. Here we present a simple set of protocols using the Targefect-HUVEC transfection kit to transfect HUVECs, HMVECs, and different types of primary endothelial cells. The kit consists of Targefect F-2, a non-lipid cationic transfection reagent with low toxicity along with two enhancer formulations- Virofect™ and a Peptide Enhancer. Virofect enhances gene transfer by using adenoviral receptors on the cell surface to enhance intracellular delivery of transfection complexes. Once the transfection complex is internalized, Virofect helps the transfection complex escape degradation in the lysosome and enhances the both duration and efficiency of transgene expression. The Peptide Enhancer increases efficiency of transgene expression by escorting genes to the nucleus. Using the protocols mentioned here researchers have been able to achieve 80-95% transfection efficiencies in HUVECs and efficiencies ranging from 40-80% in different types of primary endothelial cells (see attached figures showing data from other lab groups and a list of product citations).

### I. STORAGE

Store the transfection reagents at the following temperatures immediately upon arrival. The reagents are stable for 1 year.

Reagent	Volume	Temp.
Targefect-HUVEC (Targefect F2)	600 µl	4°C
Peptide Enhancer (PE)	600 µl	4°C
Virofect	600 µl	-20°C


**Mix each reagent prior to use by inverting 10 times.**

The Targefect-HUVEC endothelial cell transfection kit contains 3 components. Optimal transfection protocols are provided in the kit to carry out successful transfection of the following endothelial cells with high efficiency/viability.


Cell Types	Transfection Efficiency (%)
HUVECs (human umbilical vein endothelial cells)	70-95
Human dermal microvascular endothelial cells	40
Human lung microvascular endothelial cells	70-90
Human aortic endothelial cells	60
Human coronary artery endothelial cells	65
Human pulmonary artery endothelial cells	50
Human brain microvascular endothelial cells	
Bovine aortic endothelial cells	70
Bovine coronary artery endothelial cells	60
Porcine aortic endothelial cells	60
Rat aortic endothelial cells	60
Rat brain microvascular endothelial cells	

### II. PREPARATION OF CELLS

Plate Endothelial Cells at a density of 25,000-30,000 cells per cm<sup>2</sup> in the Antibiotics-Free Growth Medium. Transfect cells when cell density reaches 70-80% confluency. Preferred media: MCDB131 from VEC technologies, Media 199 plus supplements and 20% fetal calf serum or EBM from Cambrex plus supplements and 10% fetal calf serum. Cells grown in these media appear healthier and give higher transfection efficiency.



**Figure 1: Transfection of HUVECs with a GFP expression vector using the Targefect F2 reagent plus Virofect enhancer (both components of the Targefect-HUVEC kit):** 6 µg of DNA was complexed with 12 µl of Targefect-HUVEC and 25 µl of Virofect in 0.5 ml of high glucose DMEM, and then incubated at 37°C for 20 mins to form transfection complexes. 0.5 ml of transfection complexes were added to 2 ml of fresh EBM (Cambrex), 10% FCS (GIBCO) and supplements (CAMBREX) in one 60cm dish of HUVECs. The media the next day with 3 ml of complete media Transfection efficiency approx. 80%. **Data courtesy of Dr. Michael Potente, Department of Cardiology, University of Frankfurt, Germany.**



**Figure 2: Transfection of human lung microvascular endothelial cells with a green fluorescent protein expression vector using Targefect-HUVEC and the Peptide Enhancer.** Confocal images of cells transfected with a GFP-expression vector and counter-stained with rhodamine-phalloidin (actin stain). **Data**

**Courtesy of Dr. Steve Duffy and J. Murphy, UT, Southwester Medical Ctr., Dallas, TX.**

## PROTOCOL 1: (PREFERRED PROTOCOL) Transfect Endothelial Cells with Targefect-HUVEC and Peptide Enhancer

### III. FORMATION OF TRANSFECTION COCKTAIL

#### A. Preparation of DNA (Step 1)

1. Dilute the plasmid DNA with Transfection Medium serum-free high glucose DMEM or OptiMEM1) as shown in Table 1.
2. Mix thoroughly by flicking 10 times.

#### B. Preparation of transfection Complex (Step 2)

1. Mix Targefect by inverting tube 10 times.
2. Add Targefect to diluted DNA as shown in Table 1.
3. Mix thoroughly by gently flicking 10 times.
4. Mix Peptide Enhancer (PE) by inverting tube 10 times.
5. Add PE to the Targefect-DNA mixture as shown in Table 1.
6. Mix thoroughly by gently flicking 10 times.
7. Incubate the Transfection Complex at 37°C for 25 minutes.

### IV. TRANSFECTION OF CELLS

#### A. Addition of Transfection Complex (Step 3)

1. Aspirate off antibiotics-Free Growth Media from cell culture.
2. Add the appropriate amount of Transfection Complex to each well as shown in Table 2 by gently pipetting the Transfection mix along the side of the well so as not to disrupt the cells. Incubate cells with the Transfection mix in 37°C, 5% CO<sub>2</sub> humidified incubator for 2 hours.

#### B. Replacement of Transfection Complex with Antibiotic-Free Growth Medium (Step4)

1. Aspirate off the Transfection Complex from cells.
2. Add Antibiotics-Free Growth Medium to the transfected cells as shown in Table 2.
3. Incubate the transfected cells in 37°C, 5% CO<sub>2</sub> humidified incubator for 24 hours.
4. Change to Endothelial Cell Growth Medium and assay.


**Table 1: Formation of Transfection Complex**

Tissue Culture Plate	STEP 1: Preparation of DNA			STEP 2: Preparation of Transfection Complex					Total Transfection Complex (μl)
	DNA (μg)	Transfection Medium (μl) DMEM	<i>Gently Flick 10X</i>	<b>ADD Targefect (μl)</b>	<i>Flick 10X</i>	<b>ADD PE (μl)</b>	<i>Gently Flick 10X</i>	37°C for 25 mins	
96-well	0.06	60	<i>Gently Flick 10X</i>	0.25	<i>Flick 10X</i>	0.75	<i>Gently Flick 10X</i>	37°C for 25 mins	61
24-well	0.2	200		1.0		3.0			204
12-well	0.4	400		2.0		6.0			408
6-well	1.0	1000		5.0		15.0			1020

**Table 2: Transfection of Endothelial Cells with Targefect-PE Transfection Complex**

Tissue Culture Plate	Step 3: Addition of Transfection Complex			STEP 4 Replacement of Transfection Complex with Antibiotic-Free Growth Medium		
	Aspirate off Antibiotic-Free Growth Medium	<b>ADD Transfection Complex (μl)</b>	37°C 5% CO <sub>2</sub> for 1-2 hrs	Aspirate off Transfection Complex	<b>ADD Antibiotic-Free Growth Medium (μl)</b>	37°C 5% CO <sub>2</sub> for 24 hrs
96-well	Aspirate off Antibiotic-Free Growth Medium	61	37°C 5% CO <sub>2</sub> for 1-2 hrs	Aspirate off Transfection Complex	100	37°C 5% CO <sub>2</sub> for 24 hrs
24-well		204			500	
12-well		408			1000	
6-well		1020			2000	

## PROTOCOL 2: Transfect Endothelial Cells with Targefect-HUVEC and Virofect Enhancer



**Transfection complexes with Virofect are added directly to the culture media (with serum) covering the cells**

**General considerations:** Use early passage endothelial cells, avoid using collagen coated dishes as these may lower transfection efficiency. **Culture media we recommend is Media 199 with 20% serum or EBM (Cambrex) 10% FCS (Gibco) & supplements (Cambrex). We recommend using HUVEC media supplement edx with serum (10% serum concentration) if possible. If you wish to use serum-free media to culture HUVECs, then transfection complexes should be aspirated after two hours of incubation with cells.**

### Preparation of the complexes and transfection procedure:

Since reagents sometimes freeze during shipping, we recommend gently mixing the Targefect-F2 solution once upon receipt. The Targefect-F-2 reagent should be stored at 4°C. Do not vortex the Targefect-HUVEC reagent. The Peptide Enhancer can be stored at 4°C. The Virofect Enhancer should be stored at -20°C or -70°C.

**TABLE 2**

Tube #	High Glucose DMEM (Serum free)	DNA	Targefect HUVEC	Virofect Enhancer Reagent
1	1 ml	1 µg	2 µl	4 µl
2	0.6 ml	6 µg	12 µl	25 µl

Add DMEM first. Add DNA, mix well by flicking the tube about 12 times to create a vortexing action. Add Targefect next, mix well again by flicking the tube. Incubate the tubes at 37°C for 25 minutes to form the transfection complexes.

**Condition 1 above is for transfecting endothelial cells with transfection complexes in the absence of serum or any HUVEC culture media. It is very important that you try this condition as it works very well for all types of endothelial cells. At the end of the 25 minute incubation period, aspirate off the culture medium covering the cells, and add 0.2 ml of transfection complex per well of a 24-well dish. Incubate 2 hours at 37°C. Aspirate off**

**transfection complex and add 0.5 ml of complete media (with 10% serum. Incubate overnight, and replace with fresh medium the next day. Assay at 36-48 hours post transfection.**

**Condition 2 (complexes formed in tube 2)** is for transfecting cells in the presence of serum according to our fast protocol.-250 µl of transfection complex is added to 1 ml of complete media (with serum) per well of a six-well dish. The dish is swirled to enable mixing of the transfection complex with the cell culture medium and the cells are incubated at 37°C overnight and assayed for gene expression 36-48 hours post transfection.

Add 0.5 ml of transfection complex to 2 ml of complete media with serum for one 60 mm dish. For transfecting cells in a 12-well dish add 0.125 ml transfection complex to 0.5 ml of complete media.

Swirl the dish to gently mix transfection complexes with the cell culture media. Incubate overnight. Replace media the next day. Assay at 24-48 hours after transfection.

### **Citations for transfection of DNA and transfection of siRNA into endothelial cells using Targefect-HUVEC components**

#### **Citation List 1: Citations for Gene Delivery**

1. Weifei Zhu, Unni M. Chandrasekharan, Smarajit Bandyopadhyay, Sidney M. Morris, Jr., Paul E. DiCorleto, and Vikram S. Kashyap (2010)Thrombin induces endothelial arginase through AP-1 activation Am J Physiol Cell Physiol, 298: C952 - C960. Rat aortic endothelial cells
2. Bon-Hun Koo, David M. Coe, Laura J. Dixon, Robert P.T. Somerville, Courtney M. Nelson, Lauren W. Wang, Mary Elizabeth Young, Daniel J. Lindner, and Suneel S. Apte (2010) ADAMTS9 Is a Cell-Autonomously Acting, Anti-Angiogenic Metalloprotease Expressed by Microvascular Endothelial Cells Am. J. Pathol., 176: 1494 - 1504. Human Microvascular Endothelial Cells, HUVEC
3. Payal Khanna, Tara Yunkunis, Hari S. Muddana, Hsin Hsin Peng, Avery August, and Cheng Dong (2010) p38 MAP kinase is necessary for melanoma-mediated regulation of VE-cadherin disassembly Am J Physiol Cell Physiol, 298: C1140 - C1150. HUVECs
4. Catherine Chabot, Kathleen Spring, Jean-Philippe Gratton, Mounib Elchebly, and Isabelle Royal (2009) New Role for the Protein Tyrosine Phosphatase DEP-1 in Akt Activation and Endothelial Cell Survival Mol. Cell. Biol., Jan 2009; 29: 241 - 253. HUVECs
5. Arnold Johnson (2009) TNF-induced activation of pulmonary microvessel endothelial cells: a role for GSK3β. Am J Physiol Lung Cell Mol Physiol; 296: L700 - L709. Microvascular endothelial cells
6. Maud Martin, Michael Potente, Veerle Janssens, Didier Vertommen, Jean-Claude Twizere, Mark H. Rider, Jozef Goris, Stefanie Dimmeler, Richard Kettmann, and Franck Dequiedt (2008) Protein phosphatase 2A controls the activity of histone deacetylase 7 during T cell apoptosis and angiogenesis. PNAS, Mar 2008; 105: 4727 - 4732. Human umbilical vein endothelial cells (HUVECs) transfected with Targefect F-2 plus Virofect
7. Corttrell M. Kinney, Unni M. Chandrasekharan, Lori Mavrikakis, and Paul E. DiCorleto (2008) VEGF and thrombin induce MKP-1 through distinct signaling pathways: role for MKP-1 in endothelial cell migration. Am J Physiol Cell Physiol, Jan 2008; 294: C241 - C250. Human umbilical vein endothelial cells (HUVECs)
8. Primary chondrocytes
9. Smarajit Bandyopadhyay, Mohammad Z. Ashraf, Pamela Daher, Philip H. Howe, and Paul E. DiCorleto (2007) HOXA9 Participates in the Transcriptional Activation of E-Selectin in Endothelial Cells Mol. Cell. Biol; 27: 4207 - 4216. Human umbilical vein endothelial cells (HUVECs)
10. Nancy Gertzberg, Tina Gurnani, Paul Neumann, Anne-Kay Forbes, Natacha Jean-Louis, and Arnold Johnson. (2007) Tumor necrosis factor-α causes barrier dysfunction mediated by tyrosine198 and tyrosine218 in beta-actin. Am J

Physiol Lung Cell Mol Physiol, 293: L1219 - L1229.

Bovine pulmonary microvascular endothelial cells, 79% transfection efficiency using F-2 plus virofect

11. Mélanie Laramée, Catherine Chabot, Monikca Cloutier, Raphaëlle Stenne, Marina Holgado-Madruga, Albert J. Wong, and Isabelle Royal (2007) The Scaffolding Adapter Gab1 Mediates Vascular Endothelial Growth Factor Signaling and Is Required for Endothelial Cell Migration and Capillary Formation. *J. Biol. Chem.*, Mar 2007; 282: 7758 - 7769.

Endothelial cells (Targefect siRNA kit)

12. Xiaomei Qin, Jianwei Tian, Peng Zhang, Yanbo Fan, Li Chen, Youfei Guan, Yi Fu, Yi Zhu, Shu Chien, and Nanping Wang. (2007) Laminar shear stress up-regulates the expression of stearoyl-CoA desaturase-1 in vascular endothelial cells. *Cardiovasc Res*, Jun 2007; 74: 506 - 514.

Cardiovasc Res, Jun 2007; 74: 506 - 514.

Human umbilical vein endothelial cells (HUVECs)

13. Li-Wu Qian, Jianping Xie, Fengchun Ye, and Shou-Jiang Gao (2007) Kaposi's Sarcoma-Associated Herpesvirus Infection Promotes Invasion of Primary Human Umbilical Vein Endothelial Cells by Inducing Matrix Metalloproteinases. *J. Virol.*, Jul 2007; 81: 7001 - 7010.

Transfection of bacterial artificial chromosomes (BACs) into HUVECs using Targefect-HUVEC

14. Mirianas Chachisvilis, Yan-Liang Zhang, and John A. Frangos (2006) G protein-coupled receptors sense fluid shear stress in endothelial cells. *PNAS*, Oct 2006; 103: 15463 - 15468.

Bovine aortic endothelial cells , targefect F-2 plus virofect

15. Ben-Quan Shen, David Y. Lee, Hans-Peter Gerber, Bruce A. Keyt, Napoleone Ferrara, and Thomas F. Zioncheck (1998) Homologous Up-regulation of KDR/Flk-1 Receptor Expression by Vascular Endothelial Growth Factor in Vitro *J. Biol. Chem.* 273: 29979-29985.

HUVECs (Human umbilical vein endothelial cells)

16. Nanping Wang, Lynne Verna, Hai-ling Liao, Alex Ballard, Yi Zhu, and Michael B. Stemerman . Adenovirus-Mediated Overexpression of Dominant-Negative Mutant of c-Jun Prevents Intercellular Adhesion Molecule-1 Induction by LDL: A Critical Role for Activator Protein-1 in Endothelial Activation Arterioscler. *Thromb. Vasc. Biol.* 21: 1414-1420

HUVECs (Human umbilical vein endothelial cells)

17. Zhu, Hailing Liao, Nanping Wang, Kuo-Sheng Ma, Lynne K. Verna, John Y.-J. Shyy, Shu Chien, and Michael B. Stemerman LDL-Activated p38 in Endothelial Cells Is Mediated by Ras. *Arterioscler. Thromb. Vasc. Biol.* 21: 1159-1164

HUVECs (Human umbilical vein endothelial cells)

18. Nanping Wang, Lynne Verna, Neng-Guin Chen, Jasmine Chen, Hongling Li, Barry Marc Forman, and Michael B. Stemerman (2002) Constitutive Activation of Peroxisome Proliferator-activated Receptor-Suppresses Proinflammatory Adhesion Molecules in Human Vascular Endothelial Cells. *J. Biol. Chem.* 277: 34176 - 34181.

HUVECs (Human umbilical vein endothelial cells)

19. Quang-Kim Tran, D. J. Black, and Anthony Persechini (2003) Intracellular Coupling via Limiting Calmodulin *J. Biol. Chem.*, 278: 24247 - 24250.

HUVECs (Human umbilical vein endothelial cells)

20. Amir Tadros, David P Hughes, Benjamin J Dunmore, and Nicholas P J Brindle (2003) ABIN-2 protects endothelial cells from death and has a role in the anti-apoptotic effect of angiopoietin-1. *Blood*, 10.1182/blood-2003-05-1602.

HUVECs (Human umbilical vein endothelial cells)

21. Radu V. Stan.(2003) Multiple PV1 Dimers Reside in the Same Stomatal or Fenestral Diaphragm *Am J Physiol Heart Circ Physiol*, 10.1152/ajpheart.00909.2003

HUVECs (Human umbilical vein endothelial cells), HEK-293, Cos-7

22. Guifu Wu, Arjuna P. Mannam, Jiaping Wu, Simona Kirbis, Jue-Lon Shie, Christopher Chen, Roger J. Laham, Frank W. Sellke, and Jian Li (2003) Hypoxia induces myocyte-dependent COX-2 regulation in endothelial cells: role of VEGF. *Am J Physiol Heart Circ Physiol*, 285: 2420

HUVECs , HMVECs (microvascular endothelial cells)

23. Hung-Ying Kao, André Verdel, Chih-Cheng Tsai, Cynthia Simon, Henry Juguilon, and Saadi Khochbin (2001) Mechanism for Nucleocytoplasmic Shuttling of Histone Deacetylase . J. Biol. Chem., Dec 2001; 276: 47496 - 47507. HUVECs (Human umbilical vein endothelial cells)
24. Bo Lu, Yi Mu, Carolyn Cao, Fenghua Zeng, Sylke Schneider, Jiahui Tan, Jim Price, Jun Chen, Michael Freeman, and Dennis E. Hallahan (2004) Survivin As a Therapeutic Target for Radiation Sensitization in Lung Cancer. Cancer Res., 64: 2840 - 2845. HUVECs (Human umbilical vein endothelial cells)
25. Jiaping Wu, Cherie Parungo, Guifu Wu, Peter M. Kang, Roger J. Laham, Frank W. Sellke, Michael Simons, and Jian Li (2004) PR39 Inhibits Apoptosis in Hypoxic Endothelial Cells: Role of Inhibitor Apoptosis Protein-2. Circulation, ; 109: 1660 - 1667. HUVECs (Human umbilical vein endothelial cells)
26. Lingfang Zeng, Hailing Liao, Yi Liu, Tzong-Shyuan Lee, Minjia Zhu, Xian Wang, Michael B. Stemberman, Yi Zhu, and John Y-J ShyySREBP2 downregulates ABCA1 in vascular endothelial cells: a novel role of SREBP in regulating cholesterol metabolism. J. Biol. Chem., Sep 2004; 10.1074/jbc.M407817200. HUVECs (Human umbilical vein endothelial cells)
27. Radu V. Stan, Eugene Tkachenko, and Ingrid R. Niesman (2004) PV1 Is a Key Structural Component for the Formation of the Stomatal and Fenestral Diaphragms. Mol. Biol. Cell published online May 21, 2004,10.1091/mbc.E03-08-0593 HUVECs (Human umbilical vein endothelial cells), HEK-293, Cos-7 , Hela
28. Eric. T. Shinohara, Dennis Hallahan and Bo Lu (2004) The Use of Antisense Oligonucleotides in Evaluating Survivin as a Therapeutic Target for Radiation Sensitization in Lung Cancer. Biol. Proced. Online 2004;6:250-256 HUVECs (Human umbilical vein endothelial cells)
29. C. Radel and V. Rizzo (2005) Integrin mechanotransduction stimulates caveolin-1 phosphorylation and recruitment of Csk to mediate actin reorganization. Am J Physiol Heart Circ Physiol, 288: H936 - H945. HUVECs (Human umbilical vein endothelial cells)
30. Tomoaki Murakami, Hitoshi Takagi, Kiyoshi Suzuma, Izumi Suzuma, Hirokazu Ohashi, Daisuke Watanabe, Tomonari Ojima, Eri Suganami, Masafumi Kurimoto, Hideaki Kaneto, Yoshihito Honda, and Nagahisa Yoshimura (2005) Angiopoietin-1 Attenuates H<sub>2</sub>O<sub>2</sub>-induced SEK1/JNK Phosphorylation through the Phosphatidylinositol 3-Kinase/Akt Pathway in Vascular Endothelial Cells J. Biol. Chem., Sep 2005; 280: 31841 - 31849. HUVECs (Human umbilical vein endothelial cells)
31. Li-Wu Qian, Jianping Xie, Fengchun Ye, and Shou-Jiang Gao (2007) Kaposi's Sarcoma-Associated Herpesvirus Infection Promotes Invasion of Primary Human Umbilical Vein Endothelial Cells by Inducing Matrix Metalloproteinases J. Virol., Jul 2007; 81: 7001 – 7010 HUVECs (Human umbilical vein endothelial cells)
32. Smarajit Bandyopadhyay, Mohammad Z. Ashraf, Pamela Daher, Philip H. Howe, and Paul E. DiCorleto (2007) HOXA9 Participates in the Transcriptional Activation of E-Selectin in Endothelial Cells Mol. Cell. Biol., Jun 2007; 27: 4207 - 4216. HUVECs (Human umbilical vein endothelial cells)
33. Mirianas Chachisvilis, Yan-Liang Zhang, and John A. Frangos (2006) G protein-coupled receptors sense fluid shear stress in endothelial cells. PNAS, Oct 2006; 103: 15463 - 15468. HUVECs (Human umbilical vein endothelial cells)
34. Miaoliang Liu and Arie Horowitz (2006) A PDZ-binding Motif as a Critical Determinant of Rho Guanine Exchange Factor Function and Cell Phenotype. Mol. Biol. Cell, Apr 2006; 17: 1880 - 188

**List 2: References citing use of Targefect reagents (Targefect siRNA kit, catalog#0060)**

1. Boris A Hesser, Xiao Huan Liang, Gieri Camenisch, Suya Yang, David Lewin, Richard Scheller, Napoleone Ferrara, and Hans-Peter Gerber (2004) Down syndrome critical region protein1 (DSCR1), a novel VEGF target gene that regulates

- expression of inflammatory markers on activated endothelial cells. *Blood*, 10.1182/blood-2004-01-0273.(HDMVECs)  
Human Dermal microvascular endothelial cells
2. Partovian C, Simons M (2004) Regulation of protein kinase B/Akt activity and Ser(473) phosphorylation by protein kinase Calpha in endothelial cells. *Cell Signal*. 2004 Aug;16(8):951-7.
  3. Radu V. Stan, Eugene Tkachenko, and Ingrid R. Niesman (2004) PV1 Is a Key Structural Component for the Formation of the Stomatal and Fenestral Diaphragms. *Mol. Biol. Cell* published online May 21, 2004, 10.1091/mbc.E03-08-0593 HUVECs (Human umbilical vein endothelial cells), HEK-293, Cos-7, HeLa
  4. Radu V. Stan. Multiple PV1 dimers reside in the same stomatal or fenestral diaphragm. *Am J Physiol Heart Circ Physiol*, Apr 2004; 286: 1347 - 1353.
  5. Jiaping Wu, Cherie Parungo, Guifu Wu, Peter M. Kang, Roger J. Laham, Frank W. Sellke, Michael Simons, and Jian Li (2004) PR39 Inhibits Apoptosis in Hypoxic Endothelial Cells: Role of Inhibitor Apoptosis Protein-2. *Circulation*, 109: 1660 - 1667. BAECs (Bovine aortic endothelial cells)
  6. Raquel Pluvinet, Jordi Petriz, Joan Torras, Inmaculada Herrero-Fresneda, Josep M Cruzado, Josep M Grinyo, and Josep M Aran. RNAi-mediated silencing of CD40 prevents leukocyte adhesion on CD154-activated endothelial cells. *Blood*, Aug 2004; 10.1182/blood-2004-03-0817. HUVECs
  7. Seung-Hoi Koo, Hiroaki Satoh, Stephan Herzig, Chih-Hao Lee, Susan Hedrick, Rohit Kulkarni, Ronald M Evans, Jerrold Olefsky, Marc Montminy PGC-1 promotes insulin resistance in liver through PPAR-alpha-dependent induction of TRB-3. *Nature Medicine* 10, 530 - 534 (01 May 2004) Letters Hepatocytes
  8. Jeong Hoon Kim, Hongwei Li, and Michael R. Stallcup (2003) CoCoA, a Nuclear Receptor Coactivator which Acts through an N-Terminal Activation Domain of p160 Coactivators. *Molecular Cell* 2003 12: 1537-1549 MCF-7, COS-7
  9. C. Radel and V. Rizzo. Integrin mechanotransduction stimulates caveolin-1 phosphorylation and recruitment of Csk to mediate actin reorganization. *Am J Physiol Heart Circ Physiol*, Feb 2005; 288: 936 - 945. Bovine aortic endothelial cells (BAECs)
  10. Sanjeev Choudhary, Tianlin Xiao, Leoncio A. Vergara, Sanjay Srivastava, David Nees, Joram Piatigorsky, and Naseem H. Ansari. Role of Aldehyde Dehydrogenase Isozymes in the Defense of Rat Lens and Human Lens Epithelial Cells against Oxidative Stress. *Invest. Ophthalmol. Vis. Sci.*, Jan 2005; 46: 259 - 267. HLECs (Human Lens Epithelial Cells), rat and mouse lenses
  11. Yang Luo and Glenn L. Radice (2005) N-cadherin acts upstream of VE-cadherin in controlling vascular morphogenesis. *J. Cell Biol.*, Apr 2005; 169: 29 - 34. HUVEC, HDMVEC,
  12. Tomoaki Murakami, Hitoshi Takagi, Kiyoshi Suzuma, Izumi Suzuma, Hirokazu Ohashi, Daisuke Watanabe, Tomonari Ojima, Eri Suganami, Masafumi Kurimoto, Hideaki Kaneto, Yoshihito Honda, and Nagahisa Yoshimura (2005) Angiopoietin-1 Attenuates H2O2-induced SEK1/JNK Phosphorylation through the Phosphatidylinositol 3-Kinase/Akt Pathway in Vascular Endothelial Cells. *J. Biol. Chem.*, Sep 2005; 280: 31841 - 31849.
  13. Chohreh Partovian, Zhenwu Zhuang, Karen Moodie, Michelle Lin, Noriyuki Ouchi, William C. Sessa, Kenneth Walsh, and Michael Simons (2005) PKC Activates eNOS and Increases Arterial Blood Flow In Vivo. *Circ. Res.*, 97: 482 - 487
  14. Corttrell M. Kinney, Unni M. Chandrasekharan, Lori Mavrakis, and Paul E. DiCorleto (2008) VEGF and thrombin induce MKP-1 through distinct signaling pathways: role for MKP-1 in endothelial cell migration. *Am J Physiol Cell Physiol*, 294: C241 - C250.
  15. Arnold Johnson (2009) TNF-induced activation of pulmonary microvessel endothelial cells: a role for GSK3β. *Am J Physiol Lung Cell Mol Physiol*, Apr 2009; 296: L700 - L709.
  16. Catherine K. Yang, Jeong Hoon Kim, Hongwei Li, and Michael R. Stallcup (2006) Differential Use of Functional Domains by Coiled-coil Coactivator in Its Synergistic Coactivator Function with -Catenin or GRIP1. *J. Biol. Chem.*, Feb 2006; 281: 3389 - 3397.
  17. Smarajit Bandyopadhyay, Mohammad Z. Ashraf, Pamela Daher, Philip H. Howe, and Paul E. DiCorleto (2007) HOXA9 Participates in the Transcriptional Activation of E-Selectin in Endothelial Cells. *Mol. Cell. Biol.*, 27: 4207 - 4216

- 
18. Corttrell M. Kinney, Unni M. Chandrasekharan, Lin Yang, Jianzhong Shen, Michael Kinter, Michael S. McDermott, and Paul E. DiCorleto (2009) Histone H3 as a novel substrate for MAP kinase phosphatase-1 *Am J Physiol Cell Physiol*, Feb 2009; 296: C242 - C249
  19. David E. Nowak, Bing Tian, Mohammad Jamaluddin, Istvan Boldogh, Leoncio A. Vergara, Sanjeev Choudhary, and Allan R. Brasier (2008) RelA Ser276 Phosphorylation Is Required for Activation of a Subset of NF- $\kappa$ B-Dependent Genes by Recruiting Cyclin-Dependent Kinase 9/Cyclin T1 Complexes *Mol. Cell. Biol.*, 28: 3623 - 3638.
  20. Mélanie Laramée, Catherine Chabot, Monikca Cloutier, Raphaëlle Stenne, Marina Holgado-Madruga, Albert J. Wong, and Isabelle Royal (2007) The Scaffolding Adapter Gab1 Mediates Vascular Endothelial Growth Factor Signaling and Is Required for Endothelial Cell Migration and Capillary Formation *J. Biol. Chem.*, Mar 2007; 282: 7758 - 7769.
  21. Payal Khanna, Tara Yunkunis, Hari S. Muddana, Hsin Hsin Peng, Avery August, and Cheng Dong (2010) p38 MAP kinase is necessary for melanoma-mediated regulation of VE-cadherin disassembly *Am J Physiol Cell Physiol*, 298: C1140 - C1150.
  22. David Y. Lee, Jeffrey P. Northrop, Min-Hao Kuo, and Michael R. Stallcup (2006) Histone H3 Lysine 9 Methyltransferase G9a Is a Transcriptional Coactivator for Nuclear Receptors *J. Biol. Chem.*, 281: 8476 - 8485.
  23. Bon-Hun Koo, David M. Coe, Laura J. Dixon, Robert P.T. Somerville, Courtney M. Nelson, Lauren W. Wang, Mary Elizabeth Young, Daniel J. Lindner, and Suneel S. (2010) Apte ADAMTS9 Is a Cell-Autonomously Acting, Anti-Angiogenic Metalloprotease Expressed by Microvascular Endothelial Cells *Am. J. Pathol.*, Mar 2010; 176: 1494 - 1504.
  24. Weifei Zhu, Unni M. Chandrasekharan, Smarajit Bandyopadhyay, Sidney M. Morris, Jr., Paul E. DiCorleto, and Vikram S. Kashyap (2010) Thrombin induces endothelial arginase through AP-1 activation *Am J Physiol Cell Physiol*, 298: C952 - C960.
  25. Indira Padmalayam, Sumera Hasham, Uday Saxena, and Sivaram Pillarisetti (2009) Lipoic Acid Synthase (LAS1): A Novel Role in Inflammation, Mitochondrial Function, and Insulin Resistance *Diabetes*, 58: 600 - 608.
  26. Maud Martin, Michael Potente, Veerle Janssens, Didier Vertommen, Jean-Claude Twizere, Mark H. Rider, Jozef Goris, Stefanie Dimmeler, Richard Kettmann, and Franck Dequiedt (2008) Protein phosphatase 2A controls the activity of histone deacetylase 7 during T cell apoptosis and angiogenesis *PNAS*, Mar 2008; 105: 4727 - 4732.
  27. Nancy Gertzberg, Tina Gurnani, Paul Neumann, Anne-Kay Forbes, Natacha Jean-Louis, and Arnold Johnson (2007) Tumor necrosis factor- $\alpha$  causes barrier dysfunction mediated by tyrosine198 and tyrosine218 in  $\alpha$ -actin *Am J Physiol Lung Cell Mol Physiol*, Nov 2007; 293: L1219 - L1229.
  28. Catherine Chabot, Kathleen Spring, Jean-Philippe Gratton, Mounib Elchebly, and Isabelle Royal (2009) New Role for the Protein Tyrosine Phosphatase DEP-1 in Akt Activation and Endothelial Cell Survival *Mol. Cell. Biol.*, 29: 241 - 253.
  29. Young-Ho Lee and Michael R. Stallcup (2006) Interplay of Fli-1 and FLAP1 for regulation of  $\beta$ -catenin dependent transcription *Nucleic Acids Res.*, 34: 5052 - 5059.
  30. Mohini A. Patil, Susie A. Lee, Everardo Macias, Ernest T. Lam, Chuanrui Xu, Kirk D. Jones, Coral Ho, Marcelo Rodriguez-Puebla, and Xin Chen (2009) Role of Cyclin D1 as a Mediator of c-Met- and  $\beta$ -Catenin-Induced Hepatocarcinogenesis *Cancer Res.*, Jan 2009; 69: 253 - 261.



## Protocol for siRNA delivery using the Targefect-HUVEC kit

(You can also use these protocols to deliver microRNA )

The Targefect-HUVEC kit contains three components , Targefect-HUVEC (same as Targefect-F2), Peptide Enhancer and Virofect. Please find below two protocols for siRNA delivery one using the peptide enhancer and the second using the Virofect enhancer, Where ever the protocol mentions Targefect-F-2 you should use the Targefect-HUVEC component of your kit

Set up cells to be transfected so that they are about 70% confluent at the time of the experiment.

### Prepare transfection complexes as follows:

Use clear plastic tubes for complex formation. If kept frozen thaw the reagents and mix well before use (do not vortex)

Tube #	OptiMEM1	dsRNA	Targefect-HUVEC	Peptide Enhancer
1	1 ml	75 pmols	5 ul	15 ul
2	1 ml	200 pmols	5 ul	15 ul

**Make additions as follows:** Add Optimem 1 first, then add dsRNA, mix well by flicking the tube about 12 times to create a vortexing action. Add Targefect-HUVEC next, mix well again by flicking the tube and then add the Peptide enhancer, mix well again. Incubate the tubes at 37°C for 25 minutes to form the transfection complexes.

Wash cells to be transfected twice with Optimem 1. Aspirate the second wash completely. Add 1 ml of the transfection mix for 1 well for a 6-well dish (or for a 35 mm dish). Prepare 250 ul of transfection complex per well of a 12-well dish, 150 ul per well of a 24-well dish or 50 ul if using a 96-well plate Incubate the transfection complexes with the cells at 37°C for 2 hrs. Add Complete media with serum (2 ml for a 35 mm dish or one well of a 6-well plate, 1ml/well for a 12-well dish and 0.5 ml/well for a 24-well dish. Incubate overnight. Replace the media with fresh complete media the next morning and assay at 24-72 hours post-transfection.

When using fluorescently labeled dsRNA we have observed very efficient delivery into several cell types at 24 hours post transfection.

For testing for silencing of transiently transfected genes we recommend transfecting the siRNA 12 hours following transient transfection of the gene of interest and assaying for silencing of the transiently transfected gene 24 hours after transfection of the siRNA.

For testing silencing of endogenous genes we recommend analysis 48-72 hrs after transfection of the siRNA.

## Protocol for siRNA delivery using Targefect-HUVEC(Targefect-F2) plus virofect

This protocol has been standardized in HUVECs but can be applied to other cell types or suspension cells. Please note that the Targefect-HUVEC component of the Targefect-HUVEC kit is the same as Targefect-F2. Data on siRNA delivery using Targefect F-2 plus Virofect has been kindly provided by Dr Michael Potente, Department of Cardiology, University of Frankfurt, Germany

1. HUVECs used were from Cell Systems
2. The culture media is EBM with 10% FCS and supplements (Single quotes)
3. HUVEC should be 80-85% confluent at the time of transfection
4. By default we use 60 mm dishes for cell culture, therefore, the following protocol is designed for 60 mm wells
  - Add 500 ul DMEM (Optimem works as well) first
  - Add 18 ul of the desired siRNA (20 uM) and mix
  - Add 12 ul of Targefect- F2 and mix
  - Immediately add 25 ul of Virofect and mix
  - Incubate for 20 min at 37°C
  - In the meantime replace the culture media with 2 ml of fresh EBM 10% & suppl.
  - Add siRNA complexes dropwise to the cells
  - Change culture media the next day and assay for gene expression silencing after 24, 48 and 72 hours

Very recent data suggest that you can actually reduce the siRNA amount by half without losing any silencing effect. As you can see from the attached western blot, the silencing efficiency seems to greatly depend on the target sequence of the siRNA and not so much on the transfection method per se as all three methods gave comparable results. In our hands HUVECs work very well for the siRNA approach and we routinely get silencing efficiencies about 80 to 90% using different target genes. Nevertheless, the Virofect protocol is very easy to handle, suitable and time saving.

## Protocol modifications for Co-delivery of DNA and siRNA:

### Modifications suggested

To 800 ul of supernatant media covering cells in one well of a 6-well dish add 4 ul of Virofect. Incubate cells for 2 hrs. Then add 200 ul of the Transfection complex containing siRNA- Targefect F2 prepared as follows

### Preparation of transfection complex

to 200 ul of DMEM add siRNA to a concentration of 25 nM and 2 ul of Targefect F'2. Incubate 20-25 min at 37 ° C . Add the transfection complex to cells pretreated with Virofect and incubate overnight.

If transfecting DNA plus siRNA then to 200 ul of DMEM add 1 ug DNA plus siRNA (add enough to give a 25 nM concentration) . Incubate at 37 ° C for 25 min and add to cells pretreated with virofect as above.

We always recommend trying two different concentrations of siRNA - the concentration mentioned above and another which is 4 times higher.

**Note: You can also co-transfect DNA and siRNA using the Targefect-F2 reagent in combination with the peptide enhancer** If transfecting DNA plus siRNA then to prepare transfection complex mix 200 ul of DMEM add 1 ug DNA plus siRNA (add enough to give a 25 nM concentration) with 2 ul of Targefect -F2 and 4 ul of Peptide enhancer and follow instructions as in the standard protocol for siRNA delivery using the Targefect siRNA

---

**Please read the protocols accompanying the products carefully.**

**We recommend contacting tech support at 619 562 1518 for helpful tips before starting the experiment if you are a new user**

The protocols on our website need updating please do not use them. Follow the protocols that are supplied with the product

Please use antibiotic free / serum-free media for complex formation. **Regarding media for growing cells please refer to the section on “Preparation of cells” in the transfection protocol. If using complete media from Lonza or other reduced serum media for maintaining endothelial cells it is important to supplement the media with 10% serum.**

It is really important to try both protocol 1 (using the peptide enhancer) and protocol 2 (using Virofect), and to test both the conditions outlined for protocol 2 (red fonts) as in our hands this works the best. However, there does seem to be a lab to lab variation wherein one protocol works better than the other. This results partly from differences in the culture media being used. Most labs prefer the peptide enhancer protocol